Wireless Fall Detection Starter Set

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Z.o.z.

EM6795 Wireless Fall Detection Starter Set

EN

Package contents:

- EM6643 Wireless Fall Detection set (bracelet and gateway)
- EM6515 Z-wave Interface Link
- EM6505 e-Centre Base station with touch screen

Note! The RJ45 port on the gateway is NO network port. Do NOT connect the gateway to your modem/router or switch.

The potential free signal relay is not suitable to switch 230V mains power.

Connecting the devices

1. Unpack the devices
2. Connect the white data cable with RJ45 plug (already connected to the Z-wave interface link) into the alarm receiver.
3. Connect the interface link and the power adapter of the fall alarm into a mains outlet.
4. The devices are now booting up.

Using the e-Centre base station

1. Unpack the e-Centre base station
2. To get started, connect the e-Centre base station to your modem/router/switch using the provided network cable.
3. Connect the provided power adapter and connect it to the e-Centre base station, and plug it into a mains outlet.
4. The e-Centre will now boot-up.
5. After a couple of minutes you will see an empty home screen. You are now able to open the internal website of the e-Centre base station to configure it to your desire and to register it to the e-Domotica Portal to be able to receive phone calls and text messages upon the event of an alarm.
6. Open a browser on your computer or tablet and go to <http://ecentre/wiz/>. If this does not work, try to use the IP address of your e-Centre. On your e-Centre press "Devices" and go to "Settings" to find the current IP address.
7. Follow the steps on the webpage to continue

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<https://www.e-domotica.com/nl/klantenservice>

Add the interface link to your e-Centre base station

Note! If you have bought a configuration service on our website, your e-Domotica system is configured for you. You can skip this chapter.

1. Open the interface link box, using a knife or screw driver, this allows you to be able to press the Z-wave include / exclude button later.
2. Go to the 'e-Domotica' Products' page.
3. Click 'Include' to start the learn mode, the e-Centre will now listen for available Z-wave devices.
4. Press on the button of the EM6515 three times within 1.5 seconds to include it to your e-Centre.
5. The e-Centre webpage will detect an e-Domotica Sensor, a new page will appear to edit the device.
6. Enter a desired name, select a location and select an icon. Click 'Save' to continue.
7. Go to Alarm Settings to add the sensor to the personal alarm group.
8. Put the lid back on the box to close it again.

Testing and using the fall detection bracelet

The fall detection bracelet is ready to use. You can test the personal alarm system in two ways:

1. Test the fall detection alarm, drop the bracelet in front of you on the floor (height approx.. 1 meter)
 - a. After 20 seconds the bracelet will vibrate and a LED will blink. This is the first alarm period for the wearer to be able to cancel the alarm if necessary by covering the bracelet with a hand for 2 seconds.
 - b. If the alarm is not cancelled within 20 seconds the relay will switch to activate the personal alarm automatically.
2. Test the fall detection bracelet by pressing the button. The LED will blink shortly and the relay will switch to activate the personal alarm.

Disclaimer

The fall detection technology in the bracelet does not allow analysis and interpretation of all fall situations. Soft falls, slumping falls, descent controlled falls against a wall or a chair, fall down stairs etc. may not be detected by the bracelet. Considering the technology used and the target to minimize false alarms, all falls even dangerous / heavy falls may not be detected at all time, for this reason whenever the user needs assistance he/she should always press the help button on the bracelet.

Support

Find more support on our website (online manuals, instruction videos, downloads, examples etc.).

<https://www.e-domotica.com/en/customer-service>